








All models available as
 electronic Xplorer plus and
 mechanical Research plus.



Vessel Format	Type of tips	epT.I.P.S. [®]			epT.I.P.S. [®] 384	
	Pipette Electronical Mechanical	Eppendorf Xplorer [®] plus, Eppendorf Research [®] plus				
	No. of channels	4	6	8	8	12
	Volume (µL)	300 / 1,200	300 / 1,200	300 / 1,200	20 / 100	20 / 100
	Tip distance (mm)	9–33	9–20	9–14	4.5–14	4.5–9
	384 Wells (Tip distance 4.5 mm)	–	–	–	■	■
	96 Wells (Tip distance 9 mm)	■	■	■	■	■
	48 Wells (Tip distance 13 mm)	■	■	■	■	–
	24 Wells (Tip distance 19 mm)	■	■	–	–	–
	12 Wells (Tip distance 26 mm)	■	–	–	–	–
	1.5 mL / 2.0 / 5 mL Tube (Tip distance 14, 18, 20, 29 mm)	14, 18, 20, 29 mm ■	14, 18, 20 mm ■	14 mm ■	14 mm ■	–
	Agarose gel	■*	■*	■*	■	■

* Limited suitability due to volume and size of tips