

Guide pratique Ansell

Tableau des logos et symboles

Information produit

Encadré :

The diagram shows a grid of product box labels. On the left, there are four boxes: XS (pink), S (yellow), M (light pink), and XL (teal). Each box has a top section with the size letter, a middle blue section with 'POWDERED', and a bottom section with 'LATEX'. On the right, a larger purple box labeled 'L' is shown with callouts: 'Code couleur' (only on exam glove boxes), 'Taille' (the letter 'L'), 'Contenant ou non de la poudre' (the 'POWDERED' section), and 'Matière' (the 'LATEX' section).



Ce produit est composé de matériaux synthétiques et ne contient pas de latex naturel.

Données concernant la fabrication



A utiliser avant /
date de péremption



numéro de lot



Stérile, stérilisé par
radiation (rayons gamma)



Fabricant



Date de fabrication

Ansell

Données relatives au conditionnement



Une paire de gants



Environ 100 gants



Une boîte de 40 paires de gants




4 boîtes de 50 paires de gants


Instructions de conservation et manipulation




Fragile / manipuler
avec précaution



Ne pas exposer
directement au soleil



Haut - bas / cette face
vers le haut



Température limite



éviter l'exposition à l'ozone



A conserver dans un endroit
frais et sec

Informations pratiques

REF

Référence produit /
référence catalogue



Usage unique /
ne pas réutiliser



Représentant agréé au sein de
la Communauté européenne

Labels



AR 02

Gost mark (Russia)



Organisme notifié (BSI)



Kitemark
(Royaume-Uni)



Norme française