

GUIDE

Table des compatibilités de la gamme ABSolute™

Appareil de QPCR	SONDE (PROBE)				SYBR Green				
	Mix + ROX Vial	ROX Mix	Low ROX Mix	Capillary Mix	SYBR Green Mix + ROX Vial	SYBR Green ROX Mix	SYBR Green Low ROX Mix	SYBR Green Fluorescein Mix	SYBR Green Capillary Mix
ABI									
7000		•				•			
7300		•				•			
7500			•				•		
7700		•				•			
7900/7900HT		•				•			
StepOne™/StepOne™ Plus		•				•			
Bio-Rad									
iCycler™	•							•	
MyiQ™	•							•	
iQ™ 5	•				•				
Opticon™/Opticon™ 2	•				•				
Chromo 4™	•				•				
MiniOpticon	•				•				
Stratagene									
Mx4000™		•		•		•			
Mx3000P™/3005P™	•		•		•		•		
Techne									
Quantica™	•				•				
Cepheid									
SmartCycler™	•				•				
Corbett									
Rotor-Gene™	•				•				
Roche									
Lightcycler® 2.0				•					•
Lightcycler® 480	•				•				
Eppendorf									
Master Cycler Ep Realplex	•				•				
QPCR									
ABSolute Blue (Mix bleu)									
200 réactions (25 µl) 2 x 1,25 ml	-	-	-	-	-	-	-	017311	-
1600 réactions (25 µl) 16 x 1,25 ml	017285	017289	017294	-	017297	017302	017307	017312	-
400 réactions (25 µl) 1 x 5 ml	017286	017291	017295	-	017298	017303	017308	017313	-
4000 réactions (25 µl) 10 x 5 ml	-	017292	-	-	-	-	017309	017314	-
ABSolute (Mix transparent)									
200 réactions (25 µl) 2 x 1,25 ml	-	-	-	-	-	-	-	-	-
1600 réactions (25 µl) 16 x 1,25 ml	-	017527	017554	017329	-	016608	017558	017345	017349
400 réactions (25 µl) 1 x 5 ml	-	017528	017555	-	-	017337	017559	017346	-
4000 réactions (25 µl) 10 x 5 ml	-	017529	-	-	-	-	017561	-	-
QRT-PCR									
Verso 1-Step (Mix bleu)									
200 réactions (25 µl) 2 x 1,25 ml	017571	017574	017577	-	017581	017584	017587	-	-
400 réactions (25 µl) 1 x 5 ml	017573	017576	017579	-	017583	017586	017589	-	-