

Our selection ClearLine® for PCR / qPCR

- > New range of ClearLine® PCR consumables
- > Premium quality with an unbeatable price/quality ratio
- > Products in stock
- > Samples available on request

ClearLine® PCR tube strips

- > Easy to cut into individual tubes
- > Extra thin wall
- > Working volume : 0.2 ml
- > Cap strips compatible with ClearLine® PCR plates



NEW



Cat. No.	Description	Units/ case	€ excl. VAT/ case	€ excl. VAT PROMO
257052	Strip of 8 PCR tubes	250 x 8	129.30	56 -
257051	Strip of 8 optical flat caps	250 x 8	69.40	41 -
257046	Strip of 8 PCR tubes with domed cap	250 x 8	163.50	108 -
257054	8-PCR tube strips with flat cap	250 x 8	187	131 -
257053	8-PCR tube strips with flat cap	10 x 250 x 8	1164.30	913 -

ClearLine® PCR Plates

- > 96-well plate
- > Natural
- > Alphanumeric index printed in black



NEW



Cat.No	Skirt	Working Vol.	Format	Other features	Units/ case	€ excl. VAT/ case	€ excl. VAT PROMO
257047	Non skirted	200 µl	Universal	Flange around each well for easy sealing	100	223.40	143 -
257048	Skirted	100 µl	Universal	For low profile block, CFX96 compatible	100	223.40	143 -
257049	Half-skirted	200 µl	Universal / ABI	With plate flanges for alignment with purification plates	100	223.40	143 -
257050	Half-skirted	200 µl	Universal / ABI	With edges in the 4 edges of the plate	100	223.40	143 -

ThermalSeal® RTS™ ClearLine® Film

- > Applications: qPCR, long-term storage, crystallography and HPLC
- > Made of 50.8 µm polyofil
- > Operating temperature : -80 °C to +105 °C
- > For all types of plates
- > Dimensions without tab: 116.2 x 78.3 mm
- > With two detachable tabs (except roll film)
- > Optically clear, low autofluorescence
- > Good chemical resistance
- > Silicone adhesive is released only when pressure is applied for a perfect adhesion without interference with the sample



Cat.No	Description	Units/ case	€ excl. VAT/ case	€ excl. VAT PROMO
765246	ThermalSeal Film® RTS™ ClearLine®	100	301.80	195 -
765279	ThermalSeal Film® RTS™ ClearLine®	25	85.90	66 -
765265	ThermalSeal Film® RTS™ ClearLine® stérile	100	203.50	131 -
765249	ThermalSeal Film® RTS™ ClearLine® in 65 m rolls	1	997.20	717 -

Applicator

Cat. No.	€ excl. VAT/ case	€ excl. VAT PROMO
760233	54.30	22 -



Films ThermalSeal® RT2™ and RT2RR™ ClearLine®

- Application : qPCR
- Made of 50.8 µm polyester
- Operating temperature : -40 °C to +120 °C
- With two detachable tabs
- Optically clear, low autofluorescence
- High optical quality films with excellent adhesion

Cat.No	Description	Dimensions (mm)	Units/ case	€ excl. VAT/ case	€ excl. VAT PROMO
760236	ThermalSeal® RT2RR™ ClearLine® Film for non-skirted plate	118.1 x 77.8	100	256.60	164 -
760208	ThermalSeal® RT2™ ClearLine® Film for skirted or semi-skirted plate	121.9 x 79.4	100	411.10	164 -




Enter the Cat. No. on dutscher.com to view the complete Clearline® film guide

AlumaSeal® 96 ClearLine® Films


- Applications : PCR, long-term storage
- 38.1 µm aluminium
- Operating temperature : -80 °C to +120 °C
- For all types of plates
- With a tab
- Drillable
- Opaque
- Reinforced sealing for 96-well plates

Cat. No.	Description	Dimensions (mm)	Units/ case	€ excl. VAT/ case	€ excl. VAT PROMO
760231	AlumaSeal® 96 ClearLine® Film for 96 plate	117.5 x 77.8	100	97.20	62 -



ClearLine® Heat sealable films

- Available in 2 formats:
 - sheet format for manual or semi-automatic heat sealers
 - roll format for automatic heat sealers, Agilent VII Plateloc, REMP (LHS / SHS), ALPS 5000 and IntelliXseal

Photo	Model	Description	Applications	Operating temperature	Compatible plates	Cat.No	Packaging	Dimensions	€ excl. VAT	€ excl. VAT PROMO
	qPCRclear NPeel	Optical quality and peelable film	qPCR, short-term storage	-80 °C to 80 °C	Polypropylene (PP), polyethylene (PE), polystyrene (PS), cyclic olefin copolymer (COC)	257114	1 roll	610 m x 78 mm	650	500 -
						257116	100 sheets	125 x 78 mm	60	46 -
						257115	1 roll	500 m x 115 mm	636	489 -

Discover our entire range of ClearLine® adhesive and heat-sealable films at www.dutscher.com



Find our sample request form on our website in the service tab, in order to test our products for free.

Dominique DUTSCHER SAS
2C, rue de Bruxelles -
67170 BERNOLSHEIM
Tel. 03 88 59 33 90
Fax. 03 88 59 33 99
www.dutscher.com